

VTK-Unity

Medical Volume Viewer

v2.0.0 - [Google Docs Link](#)

Contact: kitware.eu/contact/

Bug tracker: [VTKUnity-BugTracker](#)

Table of contents

Limitations

Contact

Example Scene description

Volume Rendering

Image Slice

Window/Level

How to use

Monobehaviour scripts components

License



unity





ActiViz.NET

Overview

The role of this free asset is to demonstrate the [Visualization Toolkit \(VTK\)](#) rendering and processing capabilities inside a Unity application. It provides the minimal infrastructure required to create a medical image viewer in Unity. The VTK scene is rendered into Unity's rendering pipeline using Unity [low-level native plug-in interface](#) and leveraging the use of [command buffers](#).

The asset is organized with the following directory structure:

- **Plugins:** Native and managed code exposing a subset of the VTK C++ API in Unity.
- **Scripts:** [MonoBehaviour](#) scripts used to synchronize and update VTK rendering.
- **StreamingAssets:** Data files necessary for the demo scene. You might need to **copy** the contents of this inner directory to the **StreamingAssets** directory located in the **Assets** folder of your project.
- **Scenes:** Example scene using input files from the StreamingAssets directory to showcase VTK visualization and computation capabilities. An exhaustive description of the features available in the demo scene is given in section 4.



Limitations

- **OS support:** The native plugin has been implemented and tested for Windows only. Support for Linux, MacOS, and Android might come in future versions.
- **Unity Graphics API:** VTK relies on the OpenGL rendering backend, thus Unity [graphics API](#) must be set to OpenGLCore.
- **Missing features:** This asset only provides a small subset of the VTK features available in C++. While future versions of this asset plan to extend the list of available features, it is possible to access the whole set of features offered by VTK by using [Activiz](#) which provides VTK API in C#. Please contact kitware@kitware.fr for details.
<https://assetstore.unity.com/packages/essentials/tutorial-projects/vtkunity-activiz-163686>
- **Virtual and Augmented Reality:** Although VTK and ActiViz natively support rendering in VR and AR headsets such as the Hololens 2, the Unity support for OpenXR does not work with the OpenGLCore rendering backend. For this reason, the plugin only supports rendering in VR headsets using the Unity legacy VR support that is based on OpenVR and available in **Unity 2019**.

Contact

Let us know how to improve Unity support so it better fits your needs at kitware@kitware.fr, or use the [bug tracker](#) to report issues.



© 2008-2023 Kitware, Inc. & Kitware Europe
<https://www.kitware.fr/>

Example Scene description

The following features are presented in the example scene shipped with this asset and can be accessed through the VTKSceneObject.cs script. The public variables of the script are exposed in the Unity Editor to control the different features offered by the example scene. An exhaustive description is given below.

Volume Rendering

[Volume Rendering](#) is used to display a series of 2D slices as a 3D volume.



- **Volume File Name:** Specify the filename of the 3D input image (volume) to visualize. The filename can be relative to the StreamingAssets folder located in the Assets folder of the Unity project. The following extensions are supported as input:
 - VTK XML Image (.vti)
 - MetaImage (.mha/.mhd)
 - DICOM folder (Path to a folder containing DICOM images)
 - NifTI (.nii, .nii.gz)
 - Nrrd (.nrrd, .nhdr)
- **Volume Cropping:** Move the cropping sliders to control the amount of cropping applied to the volume in each anatomical direction. A value of 0.0 will disable cropping in every direction, while a value of 1.0 results in the volume being fully cropped for this specific direction.





Image Slice

It is possible to display the anatomical planes of the input volume.

The slice index can be changed interactively by moving the “Slice” slider. A value of 0.0 will compute and render the first slice in the current anatomical direction, while a value of 1.0 will present the last slice.


The current anatomical direction of the slice can be changed using the dropdown button of the “Orientation” parameter.



Window/Level

Window/Level, also known as windowing or contrast modification, defines the mapping of the input images grayscale values to black and white colors.

The brightness of the image is adjusted using the window level, while the contrast is adjusted using the window width.



How to use

To use the plugin, create a new Unity 3D project and proceed as follow:

- Change rendering backend to OpenGLCore
 - Edit->Project Settings->Player->Other Settings
 - Untick "Auto Graphics API for Windows"
 - Remove current API and add "OpenGLCore"

To try the sample scene:

- Import the VTKUnity-MedicalVolumeViewer package and open the sample scene in the Scene folder.

To start a project from scratch:

- Copy the Plugin folder into the Assets folder of your project
- Copy the Scripts folder into the Assets folder of your project
 - Add VTKCamera.cs to the main Camera of your Unity project
 - Add VTKLight.cs to every light component in your Unity project
 - Add VTKSceneObject.cs to an empty game object in the Unity scene
- Copy input images to be visualized in the StreamingAssets folder located in the Assets folder of the unity project.

Monobehaviour scripts components

The following script components are provided for users to set up a VTK scene within their Unity project.

Scripts within the Core folder are provided to synchronize the two scenes:

- *VTKCamera.cs*: This script should be attached to Unity's Camera to trigger VTK rendering. It uses Unity's [command buffers](#) to call the VTK rendering callback implemented in the native plugin. It exposes a `vtkRenderer` object in which the VTK scene objects can be rendered.
- *VTKLight.cs*: This script should be added to a Unity light object to set up the corresponding VTK scene light and synchronize light parameters. VTK supports illumination with Directional, Point and Spot light types. Shadows, indirect lighting and cookies are unsupported for now.
- *VTKSceneObject.cs*: The *VTKSceneObject.cs* script implements a `MonoBehavior` subclass that provides callback functions being called for each camera having a *VTKCamera* component. Such callbacks allow users to easily add and remove scene objects from the associated `vtkRenderer`.
- *VTKSceneCamera.cs*: This script can be attached to an empty game object to trigger rendering in the editor without starting the game. This requires all scripts to use the [ExecuteInEditMode](#) keyword.
- *VTKNativePlugin.cs*: This script provides helper functions called internally to synchronize the VTK render passes with the Unity render pipeline.
- *VTKNativePluginLite.cs*: Additional helper functions exposing a subset of the VTK C++ API in Unity in order to manage the VTK rendering stack objects such as the camera.
- *VTKNativePluginLiteScene.cs*: This script exposes a subset of the VTK C++ API in Unity in order to add, remove, and manage the VTK scene elements. The full description of features provided in the sample scene is available in the next section.

License

VTK is an open-source toolkit distributed under the OSI-approved BSD 3-clause License. See the following Copyright.txt for details.

Unset

```
/*=====
```

```
Program: Visualization Toolkit  
Module: Copyright.txt
```

```
Copyright (c) 1993-2015 Ken Martin, Will Schroeder, Bill Lorensen  
All rights reserved.
```

```
Redistribution and use in source and binary forms, with or without  
modification, are permitted provided that the following conditions are met:
```

- * Redistributions of source code must retain the above copyright notice,
this list of conditions and the following disclaimer.
- * Redistributions in binary form must reproduce the above copyright notice,
this list of conditions and the following disclaimer in the documentation
and/or other materials provided with the distribution.
- * Neither name of Ken Martin, Will Schroeder, or Bill Lorensen nor the names
of any contributors may be used to endorse or promote products derived
from this software without specific prior written permission.

```
THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS ``AS IS''  
AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE  
IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE  
ARE DISCLAIMED. IN NO EVENT SHALL THE AUTHORS OR CONTRIBUTORS BE LIABLE FOR  
ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL  
DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR  
SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER  
CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY,  
OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE  
OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.
```

```
=====*/
```